


CLINICAL SIGNIFICANCE :

Gamma -GT plays an important role in amino acid transport in the course of glutathione metabolism. The enzyme present in the serum is mainly of hepato-biliary origin. Increased enzyme activities are found in association with chronic alcoholism, different toxic liver damages, intra- and extrahepatic cholestasis, acute viral hepatitis, pancreatitis, neoplastic diseases of the liver and pancreas myocardial infarction as well as with diabetes mellitus.

TEST PRINCIPLE :

Optimized kinetic determination of g-glutamyltransferase (g-GT). The increase of the absorbance at 405 nm, due to the formation of the 5-Amino-2- nitrobenzoic acid, is proportional to the g-GT-activity.



REAGENTS COMPOSITION :

R1: Glycylglycine 150 mmol/l

R2: L-Gamma-glutamyl-3- carboxy-4-nitroanilide 6 mmol/l

KIT CONTENTS :

	CODE No. GGT02 (5x10 ml)	CODE No. GGT03 (5 X 25 ml)
Pack size :		
Reagent 1	5 x 8 ml	5 x 20 ml
Reagent 2	1 x 10 ml	1 x 25 ml

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2 - 8 °C, and contaminations prevented during their use. Do not use reagents over the expiration date. Once opened the reagent is stable for 1 month On-board the analyser at approximately 10°C.

SAMPLES : Serum or heparinised Plasma, Fasting preferable
 Serum should be separated from blood as soon as possible.
 Sample must be free of hemolysis

INTERFERENCES

The following analytes were tested up to the levels indicated and found not to interfere with Ascorbic Acid up to 30 mg/dL, Bilirubin up to 40 mg/dL, Hemoglobin up to 500 mg/dL and Triglycerides up to 2000 mg/dL.

ASSAY CONDITIONS:

Wavelength :	405 nm
Cuvette:	1 cm light path
Constant temperature	37°C
Reaction (Mode).....	Kinetic
Kinetic Factor.....	2608
Delay	60 sec
Read time.....	180 sec
Linearity.....	1800 U/L
Unit.....	U/L
Blanking.....	D. Water
Slope of reaction.....	Increasing


PROCEDURE :

Pipette into test tubes labeled Test (T) as follows:

	T
Reagent R1	400 l
Reagent R2	100 l
Mix and wait for 5 mint	
Specimen	25 l

Mix and read the change of absorption (A) with delay time 60 sec and read time 180 sec at 37°C

Assay Procedure summary:



CALCULATIONS :

GGT Activity in U/L = Abs of Test x 2608

REFERENCE RANGE :

	Female	Male
Adults:	9 - 36	12 - 64
Children:		
1 Day– 6 Months	15 - 132	12 - 122
6 Months – 1 Year	1 - 39	1 - 39
1 – 12 Years	4 - 22	3 - 22
13 – 18 Years	4 - 24	2 - 42

The above reference range is guideline and all the laboratories must establish their own Age sepcific normal reference range. Final diagnosis should be made with correlation of clinical factors.

PRECAUTIONS :

1. Storage conditions as mentioned on the kit to be adhered.
2. Use clean glassware and microtips while pipetting
3. Avoid contamination of the reagent during the assay process.
4. Before the assay begins, bring all the reagents to room temperature.
5. If a larger volume of reagent is required for the absorbance reading, requisite volume can be taken in multiples, keeping the same ratio of reagent to specimen
6. Do not freeze or expose the reagents to high temperature and protect from direct sunlight as it will affect the performance of the kit.
7. Programmes for specific autoanalysers are available on request.
8. For accuracy of results, the assay procedure, reagent preparation and storage has to be meticulously followed.
9. As with all the diagnostic procedures, the physician should evaluate data obtained by the use of this kit in light of other clinical information.

LINEARITY AND DETECTION LIMIT :

The assay is linear up to GGT activity upto 1800 IU/L. The results of the performance characteristics depend on the analyzer used. If the results obtained were greater than linearity limit, dilute the sample 1 : 5 with Normal Saline and multiply the result by 5.

BIBLIOGRAPHY

1. Szasz, G. Kinetic photometric method for serum gamma glutamyltranspeptidase. Clin. Chem., 15 (1969) 124.
2. Shaw LM, Stromme JH, London JL, Theodorsen L. International Federation of Clinical Chemistry, (IFCC), Scientific Committee, Analytical Section. IFCC methods for the measurement of catalytic concentration of enzymes. Part4. IFCC method for gamma-glutamyltransferase [(gamma-glutamyl)-peptide: amino acid gamma-glutamyltransferase, EC 2.3.2.2]. J Clin Chem Clin Biochem 1983;21:633-46.
3. Thomas L. Clinical Laboratory Diagnostics. 1st ed.Frankfurt: TH-Books Verlagsgesellschaft;1998.p.80-6.