


**CLINICAL SIGNIFICANCE :**

Gamma -GT plays an important role in amino acid transport in the course of glutathione metabolism. The enzyme present in the serum is mainly of hepato-biliary origin. Increased enzyme activities are found in association with chronic alcoholism, different toxic liver damages, intra- and extrahepatic cholestasis, acute viral hepatitis, pancreatitis, neoplastic diseases of the liver and pancreas myocardial infarction as well as with diabetes mellitus.

**TEST PRINCIPLE :**

Optimized kinetic determination of g-glutamyltransferase (g-GT). The increase of the absorbance at 405 nm, due to the formation of the 5-Amino-2- nitrobenzoic acid, is proportional to the g-GT-activity.



**REAGENTS COMPOSITION :**

R1: Glycylglycine 150 mmol/l

R2: L-Gamma-glutamyl-3- carboxy-4-nitroanilide 6 mmol/l

**KIT CONTENTS :**

CODE No. : GGTLW01  
 Pack size : (5 x 10 ml)  
 Reagent 1 4 x 10 ml  
 Reagent 2 2 x 5 ml

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2 - 8 °C, and contaminations prevented during their use. Do not use reagents over the expiration date. Once opened the reagent is stable for 1 month On-board the analyser at approximately 10°C.

**SAMPLES :** Serum or heparinised Plasma, Fasting preferable  
 Serum should be separated from blood as soon as possible.  
 Sample must be free of hemolysis

**Materials Required but not provided :**

Gamma GT IFCC Calibrator and Control  
 (Use of assayed QC sera is recommended to validate test result).

**INTERFERENCES**

The following analytes were tested up to the levels indicated and found not to interfere with Ascorbic Acid up to 30 mg/dL, Bilirubin up to 40 mg/dL, Hemoglobin up to 500 mg/dL and Triglycerides up to 2000 mg/dL.


**ASSAY CONDITIONS:**

Wavelength : ..... 405 nm  
 Cuvette: ..... 1 cm light path  
 Constant temperature ..... 37°C  
 Reaction (Mode)..... Kinetic  
 Calibrator Conc..... See calibrator value  
 Delay ..... 60 sec  
 Read time..... 180 sec  
 Linearity..... 800 U/L  
 Unit..... U/L  
 Blanking..... D. Water  
 Slope of reaction..... Increasing

**PROCEDURE :**

- Mix 13 µl Sample / Calibrator with 200 µl R1, mix and wait for 5 mints.
- Mix 50 µl R2 and read initial absorbance after 60 sec. and measure final absorbance after 180 sec against distilled water blank at 405 nm by an Analyser.
- Calculate absorbance change per minute ( $\Delta A/\text{min}$ ).

**Assay Procedure summary:**



**CALCULATIONS (With Calibrator):**

$$\text{GGT (U/L)} = \frac{\Delta \text{Abs} / \text{mint of sample}}{\Delta \text{Abs} / \text{mint of Calib.}} \times \text{Conc. of Calib.}$$

**CALCULATIONS (With Factor at 37°C):**

$$\text{GGT Activity in U/L} = \Delta \text{Abs/mint of Test} \times \text{KF}$$

**Gamma GT (GGT) Calibrator (IFCC) is recommended to use for calculation of Kinetic Factor in LW SERIES and other fully automated Analysers. Kinetic Factor (KF) value is Lot Specific.**

**REFERENCE RANGE :**

	Female	Male
Adults:	9 - 36	12 - 64
Children:		
1 Day– 6 Months	15 - 132	12 - 122
6 Months – 1 Year	1 - 39	1 - 39
1 – 12 Years	4 - 22	3 - 22
13 – 18 Years	4 - 24	2 - 42

*The above reference range is guideline and all the laboratories must establish their own Age sepcific normal reference range. Final diagnosis should be made with correlation of clinical factors.*

**PRECAUTIONS :**

- Storage conditions as mentioned on the kit to be adhered.
- Avoid contamination of the reagent during the assay process.
- If a larger volume of reagent is required for the absorbance reading, requisite volume can be taken in multiples, keeping the same ratio of reagent to specimen
- Do not freeze or expose the reagents to high temperature and protect from direct sunlight as it will affect the performance of the kit.
- Programmes for specific LW autoanalysers are available on request.
- For accuracy of results, the assay procedure, reagent preparation and storage has to be meticulously followed.
- As with all the diagnostic procedures, the physician should evaluate data obtained by the use of this kit in light of other clinical information.

**LINEARITY AND DETECTION LIMIT :**

The assay is linear up to GGT activity upto 800 IU/L. The results of the performance characteristics depend on the analyzer used. If the results obtained were greater than linearity limit, dilute the sample 1 : 5 with Normal Saline and multiply the result by 5.

**BIBLIOGRAPHY**

- Szasz, G. *Kinetic photometric method for serum gamma glutamyltranspeptidase. Clin. Chem., 15 (1969) 124.*
- Shaw LM, Stromme JH, London JL, Theodorsen L. *International Federation of Clinical Chemistry, (IFCC), Scientific Committee, Analytical Section. IFCC methods for the measurement of catalytic concentration of enzymes. Part4. IFCC method for gamma-glutamyltransferase [(gamma-glutamyl)-peptide: amino acid gamma-glutamyltransferase, EC 2.3.2.2]. J Clin Chem Clin Biochem 1983;21:633-46.*
- Thomas L. *Clinical Laboratory Diagnostics. 1st ed.Frankfurt: TH-Books Verlagsgesellschaft;1998.p.80-6.*