

**MICROALBUMIN (Albumin in urine)**

In vitro diagnostic reagents for the quantitative determination of albumin in urine (MAU) by means of particle-enhanced turbidimetric immunoassay.

**Clinical Significance**

Increased albumin excretion detectable only by sensitive immunoassay (microalbuminuria) has been used for some years as a predictor of incipient nephropathy and cardiovascular disease in diabetic patients. Microalbuminuria has also been associated with hypertension and increased risk of cardiovascular disease in non-diabetic patients. Microalbuminuria occurs in response to acute inflammatory conditions such as ischaemia, trauma, and thermal injury, surgery, pancreatitis, and inflammatory bowel disease. In many of these conditions albumin excretion increases within minutes or hours of the initiating stimulus and only last for 24 to 72 h. The degree of microalbuminuria is proportional to the severity of the inflammatory insult, is predictive of outcome, and is not associated with any other features of renal impairment. Interest in measuring subclinical elevations in the albumin excretion rate has focused on individuals with an already established diagnosis of diabetes or essential hypertension. Providing proper care is taken to minimise the influence of exercise and poor metabolic control of the albumin excretion rate, the urinary albumin level has proved to be an excellent predictor of the progression to overt nephropathy in both insulin-dependent and non-insulin dependent diabetes.

**Principle**

This MAU test is based upon the reactions between albumin and latex-covalently bound antibodies against human albumin. MAU values are determined photometrically.

**Reagents**

**A - Buffer** - 45 mL of phosphate buffer, pH: 8.5, < 0,1 % sodium azide as preservative.

**B - Latex reagent** - 5 mL of a suspension of latex microparticles covalently bound anti-albumin antibodies in a neutral aqueous solution, and < 0,1 % sodium azide as preservative.

**C - Calibrator** - 1 mL. Human - based reference fluid. Preservative: sodium azide, 0.075 %. All raw materials of human origin used in the manufacture of this product showed no reactivity when tested for HBsAg, anti-HIV-1/2 and HCV with commercially available test methods. However, this product should be handled as though capable of transmitting infectious diseases

**Reagent Preparation**

Working Reagent is prepared with 1 part of Latex Reagent and 9 parts of Buffer Reagent. Prepare a fresh WR based on its workload. Shake gently the reagents before pipetting.

**Calibration Curve and Controls**

Analytical Range up to 250 mg/L.

Calibrator 1	100 µl of MAU Calibrator*
Calibrator 2	100 µl of Calibrator 1 + 100 µl of Saline Solution
Calibrator 3	100 µl of Calibrator 2 + 100 µl of Saline Solution
Calibrator 4	100 µl of Calibrator 3 + 100 µl of Saline Solution
Calibrator 5	100 µl of Saline Solution

(\* ) See values on the label or on the insert. Multiply by the appropriate factor.

For quality control use a suitable control material. The control intervals and limits must be adapted to the individual laboratory requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits. Control must be assayed and evaluated as for patient samples.

**Stability**

Reagents in the original vial is stable to the expiration date on the vial label when capped and stored at +2 - +8°C. Immediately following the completion of an assay run, the reagent vial should be capped until next use in order to maximize curve stability. Once opened the reagent can be used within 1 month if stored tightly closed at +2 - +8°C after use. Do not freeze reagents. The MAU latex reagent should have a white, turbid appearance free of granular particulate. Visible agglutination or precipitation may be a sign of deterioration, and the reagent should be discarded. The MAU buffer reagent should be clear and colourless. Any turbidity may be sign of deterioration and reagent should be discarded. WR is stable for up to two weeks at 4°C. **It is recommended that each Laboratory prepares a fresh Working Reagent based on its workload.**

**Precautions**

For in vitro diagnostic use only. Do not pipette by mouth. Reagents containing sodium azide must be handled with precaution. Sodium azide can form explosive azides with lead and copper plumbing. Since absence of infectious agents cannot be proven, all specimens and reagents obtained from human blood should always be handled with precaution using established good laboratory practices. Disposal of all waste material should be in accordance with local guidelines. As with other diagnostic tests, results should be interpreted considering all other test results and the clinical situation of the patient.

**Materials required but not provided**

Spectrophotometric analyser. Controls.

**Specimens**

Use 12 or 24 hour collection. centrifuge urine specimens. Screen these specimens using an albumin test strip. If the result is negative (approx. below 300 mg/L), analyse the specimens undiluted. If the result is positive, dilute the specimen with specific protein sample diluent to obtain a concentration below 250 mg/L .

**Procedure**

Wavelength	600 nm		
Temperature	37°C		
Cuvette	1cm light path		
Measurement against distilled water blank.			
Bring the reagents at 37°C and pipette:			
	Calibrator	Sample	Blank
Calibrator	3 µl	---	---
Sample	---	3 µl	---
Distilled Water	---	---	3 µl
Work. Reagent	500 µl	500 µl	500 µl
Mix and measure absorbance immediately (A1) incubate 4 min (37°C), after incubation read absorbance (A2).			

**Calculation**

Plot the calibration curve and the sample concentration is obtained by interpolation the sample absorbance (A2-A1) in the calibration curve. If is an one point calibration:

$$\frac{(A2 - A1)_{\text{sample}} - (A2 - A1)_{\text{blank}}}{(A2 - A1)_{\text{calibrator}} - (A2 - A1)_{\text{blank}}} \times \text{Calibrator Concentration}$$

**Linearity**

The range interval for the multipoint calibration method is from 0 to 250 mg/L. With this method you can use the one point calibration procedure using a calibrator without dilutions, because it is linear at least up to 170 mg/L. When values exceed the range the samples should be diluted with saline solution and the result should be multiplied by the appropriated factor. In an one point calibration, when values exceed 125 mg/L, the samples should be diluted with saline solution and the result should be multiplied by the appropriated factor.

**Reference Values**

For timed overnight urine collections an albumin excretion rate greater than 20 µg/min is considered to abnormal. These data are to be interpreted as a guide. Each laboratory should establish its own reference intervals.

**Sensitivity**

Calculating the mean plus 3SD of twenty replicates of zero standard resulted in a lower limit of detection less than 5 mg/L.

**Prozone Effect**

The system did not show prozone phenomenon at least up to 400 mg/L.

**Assay Precision**

Intra-assay coefficients of variation (CV) for three samples (MAU values ranging from 30 to 150 mg/L) were between 1.9 and 3.7 %. Inter-assay CVs were between 2.4 and 4.5 %.

**Method comparison**

70 samples were correlated with a commercial procedure. When comparing the results by lineal regression the result was:  $y = 0,96 x + 4,1$  and  $r = 0,981$

**Analytical characteristics have been obtained in a single experiment in a conventional spectrophotometer. Therefore, the data expressed in the present document should be interpreted as a guide example.**

**Literature**

Winocour PH. Microalbuminuria Brni 1992;1 304;1196-7 Marshall SM. Screening for microalbuminuria: which measurement. Diabetic Medicine 1991; 8: 706-11  
 Osherg Y et al. Effects of storage time and temperature on measurement of small concentrations of albumin in urine. Clin Chem 1990; 36:1428-30  
 Gosling P. Microalbuminuria: a sensitive indicator of non-renal disease?. Ann Clin Biochem 1995; 31439-41  
 Passing H, Bablok W. A new biometrical procedure for testing the equality of measurements from two analytical methods. Application of linear regression procedures for method comparison studies. Part I. J Clin Chem Clin Biochem 1983; 21:709-20.  
 Sonderdruck aus DG Klinische Chemie Mitteilungen 1995; 26: 207 - 224

Significados de los símbolos indicados en las etiquetas. Explanation of symbols used on labelling. Explication des symboles figurant sur les étiquettes. Spiegung der symbole utilizzati sull'etichetta. Significado dos símbolos indicados nas etiquetas. Erläuterung der symbole auf den etiketten.

												
Fecha de Caducidad Expiry Date Date de Péremption Data di Scadenza Data Expiração Verwendbar bis	Temperatura de Almacenamiento Storage Temperature Temperatura de Conservazione Temperatura de Conservação Temperatura de Lagerung	Número de Lote Lot Number Número de Lote Numero di Lote Número de Lote Chargen-Nr.	Para Diagnóstico In Vitro For In Vitro Diagnostic Usage In Vitro Per Uso Diagnostico In Vitro Utilizar em Diagnostico In Vitro In Vitro Diagnostikum	Número de catálogo Catalog Number Número de catálogo Numero di catalogo Número de catálogo Katalognummer	Conformidad Europea European Conformity Conformité aux normes européennes Conformità europea Conformidade com as normas europeias CE-Konformität/Bezeichnung	Fabricado por Manufactured by Fabricado por Fabbricato da Fabricado por Hergestellt	Reactivo Reagent Réactif Reagent Reagente Reagens	Calibrator Calibrator Calibratore Calibratore Calibrador Calibrador	Buffer Buffer Tampone Tampone Buffer Puffer	LYOPH Lyophilized Lyophilisé Lyophilise Lyophilisat	Conc. Concentration Concentration Concentração Konzentration	Control H / Control L Control Alto / Control Bajo Control High / Control Low Contrôle élevé / Contrôle Bas Controlo Alto / Controlo Basso Controlo Alto / Controlo Baixo Kontrolle Hoch / Kontrolle Niedrig

