



Product features

Mini Fully Auto Chemistry Analyzer
Discrete/random access
STAT sample priority
Easy to use
User friendly interface
No computer
Internal thermal printer
Direct print report
All-in-one machine

Working Conditions

Temperature: 10°C ~ 30°C
Humidity: < 85%
Power Supply: 100~240V 50/60HZ
Water Consumption: < 5L



Reaction/Reagent/Sample Handling:

81 cuvettes
24+5(Tube Ø 12*100) sample positions
1 calibration position / 1 QC position / 1 ST position
54 positions for reagent
27 bottles with volume of 8ml
27 bottles with volume of 24ml
LED light source with 8 channel transmissometry, Rear Spectrophotometry
Extra 1 channel scattering (Option : SL120 C)
Liquid level detection / Collision protection
High pressure water flush interior the probe

Other

Test speed : Constant 120T/H
Dimension : 305mm(L) x 605mm(W) x 330mm(H)
Net weight : 22KG

 **DIATEK** HEALTHCARE PVT. LTD.
(An ISO 9001 : 2008 Company)

Corporate Office :

P-27, Kalindi Housing Scheme,
Kalindi, Kolkata, WB, India, Pin - 700089
Email : diatek@rediffmail.com,
TEL : +91-33-25223075/3448/3646
Fax : +91-33-2522401

For enquiries please
Call us Toll Free on

1800 120 5020



Dia-CHEM 120

Automatic Chemistry Analyzer



Internal printer
One button to test
Touch screen
Internal operating system

SEE it

Dia-CHEM
120

Smart Economic Efficient Automatic Chemistry Analyzer