



LW W12 Microplate Washer

- 12-way wash head interchangeable
- Programmable up to 64 wash protocols
- Minimal residual volume
- Low well-to-well variation
- Automatic self-diagnostic program
- Automatic liquid monitoring system
- Touch screen , LCD display

LANDWIND
MEDICAL

• Specifications



- **WASH HEAD**
12 heads
- **WASH MODE**
Row and plate
- **WASH ROWS**
1-12 rows
- **PRIME VOLUME**
50-700 μ l
- **MOVING CYCLE**
1-9 cycles
- **CLEANING TIME**
1-255s
- **CLEAR BOTTOM**
Yes/No
- **WASH BOTTOM**
Yes/No
- **BOTTLE CAPACITY**
>2L
- **WASH PROGRAMS**
Up to 64
- **TYPE OF PLATE**
Up to 10
- **RESIDUAL VOLUME PER WELL**
 $\leq 2 \mu$ l for V/ U bottom plates, $\leq 3 \mu$ l for flat bottom plates
- **DISPLAY**
128 \times 64 pixels with backlight
- **KEYPAD**
Touch screen
- **OPERATING ENVIRONMENT**
Input power: <80W
Power supply: AC220V/110V, 50/60Hz
Ambient temperature: 10 $^{\circ}$ C-30 $^{\circ}$ C
Relative humidity: $\leq 20\%$ -85%
Weight: 7.5Kg (Bottles not included)
Dimension: 380mm(H) \times 350mm(W) \times 190mm(D)

Shenzhen Landwind Industry Co., Ltd

Headquarter:

Landwind Science & Technology Park, Tangtou Avenue,

Shiyan Town, Bao'an District, Shenzhen China 518108

Tel: +86-755-83933788 Fax: +86-755-27353240

Email: Sales@landwindmedical.com

Web: www.landwindmedical.com

DISTRIBUTOR:

Specifications and appearance are subject to change without prior notice.

P/N: